




Available Mattresses

Ascent and Soft touch Mattresses are Available in Bamboo Cloth and Cleanable Moisture Resistant Stretch Vinyl, Offer Superior Pressure Redistribution and Support, and have a firm edge to help with transfers in and out of bed.

Available Sizes

| Twin | Full | Queen | Dual King |
|-------|-------|-------|--|
| 38x80 | 53x80 | 60x80 | 2 Twins 2x38x80 independent functions |



ASCENT 6" MATTRESS: Most supportive mattress offering zoned support and air channels to help dissipate moisture and heat

SOFT TOUCH 7" MATTRESS: Memory Foam mattress that offers support while contouring with the body



“ My new Supernal bed frame and mattress from Transfer Master has been a wonderful improvement in the quality of my life ”

Ted, SC, USA

Mobility Plus®

Transfer Master™



Supernal[®]

SLEEP SYSTEM



**MAXIMIZE YOUR
INDEPENDENCE,
TAKE CONTROL
OF YOUR COMFORT**



Features

- **10 INCH HEIGHT RANGE:** 10.5" - 20.5"
(Deck Height Without Mattress)
- **HEAD AND FOOT ELEVATION** - 65° Head Elevation, 35° Foot Elevation, and Hi-Low Adjustable
- **WIRELESS REMOTE** - Backlit wireless hand control with convenient built-in flashlight
- **EUROPEAN STYLE HEAD TILT** - Cradles your pillow w/ a European style head tilt
- **WALL-HUGGER** - Bed stays close to wall when adjusted, so you are always within reach of your nightstand
- **VARIABLE MASSAGE SETTINGS**
- **EMERGENCY POWER DOWN** - In case of power outage, the bed may be brought down to its fully flat and lowest position. This is done with the push of a button located at the foot of the bed.
- **WEIGHT CAPACITY:**
Twin: 400 Lbs **Full/Queen:** 500 Lbs

Available Accessories



Half Rails
(head section only)



Bamboo
Rail Covers



5" Locking
Caster Wheels



Battery Backup



Headboard