

There is Power Mobility. And then there's

# POWERFUL MOBILITY



Vita Monster



There is Power Mobility. And then there's

# POWERFUL MOBILITY

## Vita Monster

The Vita Monster has a double A-Arm front suspension and 3D rear suspension for a smooth ride and maximum shock absorption. It comes with a full LED light package, rearview mirrors, and a roll bar. The Vita Monster's top speed is 11.5 mph, with a range of up to 28 miles per charge. Other features for the Vita Monster include a sporty 20" captain's seat, cup holder, and adjustable tiller.



MAX SPEED  
**11.5**  
MPH



LOAD CAPACITY  
**350**  
LBS



BATTERY  
**79 Ah**  
SLA



RANGE  
**28**  
MILES

Overall Length	64" in / 162 cm
Overall Width	31" in / 79 cm
Overall Height (at Maximum Tiller Height)	54" in / 136 cm
Overall Height - at Minimum Tiller Height	46" in / 117 cm
Legroom - Floorboard Length x Width	13" in / 132 cm
Deck Height - Floorboard to Floor	9.4" in / 24 cm
Ground Clearance	4" in / 10 cm
Turning Radius	53" in / 135 cm
Weight with Std.	330 lbs / 150 kgs
Weight without Batteries	220 lbs / 100 kgs
Seat Type	Captain Seat

Tire Type	Air Tire
Rim Type Rear	Aluminum
Rim Type Front	Aluminum
OEM Spec Rear	N/A
Outside Diameter Rear	15" in / 38 cm x 6.3" in / 16 cm
Tire Width Rear	6.3" in / 16 cm
OEM Spec Front	N/A
Outside Diameter Front	13" in / 33 cm x 4.5" in / 11.5 cm
Tire Width Front	4.5" in / 11.5 cm
Battery Type & Power	12V 75Ah x 2pcs
Brake System	Electromagnetic
Climbing Angle - Maximum Recommended	12 degree

Charger Type	8 Amp, Off Board 120/240 Volt, 50/60Hz
Charging Time	8-10hrs
Controller	S-DRIVE 200 Amp
*Driving Range	28 miles / 45 km
Load Capacity	350 lbs / 160 kgs
Motor & Drive Location	900 watts
*Speed - Maximum	11.5 mph / 18.5 kph
Colors	Midnight Black
Tiller Type	D Type Handle
Throttle Control Type	Pot Meter
Anti-Tipper Wheels-Rear	Yes