

<b>Model</b>	Go Chair®
<b>Part Number</b>	GoChair
<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b> <sup>1</sup>	Up to 3.8 mph
<b>Ground Clearance</b> <sup>3</sup>	1.75" at motor
<b>Turning Radius</b> <sup>3</sup>	25.5"
<b>Overall Length</b> <sup>3</sup>	28.5" (33" with footplate)
<b>Overall Width</b> <sup>3</sup>	25.5"
<b>Front Wheels</b>	6" solid casters
<b>Drive Wheels</b>	9" solid
<b>Rear Casters</b>	3" solid anti-tips
<b>Suspension</b>	Limited
<b>Seat-to-Ground Height</b>	22" - 23"
<b>Seat-to-Platform Height</b>	16.25" - 17.25"
<b>Range Per Charge</b> <sup>1,2</sup>	Up to 13.2 miles (200 lbs.) Up to 8.7 miles (300 lbs.)
<b>Heaviest Piece When Disassembled</b> <sup>3</sup>	36 lbs. (rear section)
<b>Standard Seating</b>	Type: Compact, foldable Material: Black vinyl Dimensions: Width: 18", Depth: 17", Weight: 30 lbs.; Width: 20", Depth: 18", Weight: 32 lbs.;
<b>Dual Braking System</b>	Regenerative and electromechanical
<b>Standard Electronics</b>	40 amp, PG GC3
<b>Battery Requirements</b> <sup>4</sup>	(Two) 18 AH, 12 V, deep-cycle Weight: 14 lbs. each
<b>Battery Charger</b>	Off-board, 2 amp
<b>Warranty</b>	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Batteries: 13-month limited

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

**HCPCS Code: K0899**

**POWER MOBILITY DEVICE, NOT CODED BY DME PDAC  
OR DOES NOT MEET CRITERIA**

1. Pride<sup>®</sup> FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

0004187-032822

**GOCHAIR**<sup>®</sup>

**GET  
OUT  
THERE!**



**LIVE YOUR BEST<sup>®</sup>**

pridemobility.com |    

# GET OUT THERE

## WITH THE #1 PORTABLE POWER CHAIR!

Convenient storage pocket

Five-point swivel seat with latch release

One-hand feather-touch disassembly

FDA Class II Medical Device<sup>1</sup>



Height and width adjustable armrests

Dual swing-away, under seat storage bins

Large footplate for added comfort

The Go Chair® is #1 for portability! With convenient feather-touch disassembly, the Go Chair allows you to experience lightweight travel and independence on the go. Enjoy dual swing-away storage bins and a compact turning radius for outstanding maneuverability. Available in an array of contemporary colors, it's time to get out there with the Go Chair!

### Additional Features

- 18" x 17" or 20" x 18" standard 5-point swivel seat with latch release
- Compact turning radius
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- Large footplate
- One hand feather-touch disassembly
- 60" lap belt
- The heaviest piece is only 36 lbs.

### Available Colors



### Accessories\*

- Armrest mount cane clip
- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Lap belt (50", 70")
- Oxygen tank holder
- Rear mount bracket
- Saddlebag
- Safety flag
- Weather cover
- Inline right and left swing-away joystick
- XLR USB charger
- Additional battery pack

### Military Patches



\* Only one rear accessory can be used at a time.



Disassemble quickly for transport or storage.

Store items close at hand with convenient under seat storage.

Personalize your Go Chair® with one or more optional accessories.