



PRIDE[®]
POWER
RECLINER
COLLECTION¹

“

I LOVE MY PRIDE LIFT CHAIR!
I CAN FULLY RECLINE AND BE
COMPLETELY COMFORTABLE!
AND, THE FABRIC MATCHES
MY LIVING ROOM PERFECTLY!

– GEORGE R.



SCOTT MEUSER
CHAIRMAN & CEO
PRIDE MOBILITY®

With over three decades in the power lift recliner industry, we know that our business is about you, our customer. We're humbled when you choose to bring one of our power lift recliners into your home, knowing that it must fit your lifestyle and decor.

Pride® Power Lift Recliners is the culmination of your direct feedback. You wished for luxurious comfort, exceptional quality, state-of-the-art features and elegant designs. Over a 30-year period, we listened and refined, then listened some more. The result is the Pride line of power lift recliners designed for you, your home and your lifestyle.





We understand that bringing a power lift recliner into your home isn't just a purchase, but an investment – toward your comfort, home décor, and lifestyle. We deeply value that you, our customer, have enabled us to meet those vital goals through Pride Power Lift Recliners.

A handwritten signature in black ink that reads "Scott Meuser". The signature is fluid and cursive, with a long horizontal line extending from the end of the name.

RECLINE POSITION



FUNCTIONS AND FEATURES

-  1. The full recline position is ideal for napping or relaxing.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. Dual side pockets are ideal for storing remotes, books or magazines.
-  4. The footrest extension allows you to relax and stretch out completely.




FOOTREST EXTENSION⁵



BACKREST⁵



USB REMOTE



5. The USB remote adjusts all functions to help you find the perfect position.



6. The backrest moves up and down independently to help you find the perfect position.



ELEGANT COMFORT AND RELAXATION

STANDARD FEATURES



7. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.

6.

INFINITE POSITION⁵



7.

“

MY PRIDE RECLINER GIVES ME
THE SUPPORT I NEED
AND SO MUCH MORE!
I LOVE HOW EASY IT IS
TO ADJUST MY POSITION!

– EMILY R.

HIGH-RESILIENCY FIBER FILL | ■
for maximum comfort

FABRIC | ■
Soft, comfortable, luxurious fabrics are durable as they are stylish. Tested to commercial standards, guaranteed to last.



..... ■ | **STRONG WOOD FRAME**
Premium furniture-grade hardwood is robust, durable and long-lasting.




UNIQUE SPRING DESIGN | ■
Interlocking springs ensure maximum support while delivering individualized comfort.

..... ■ | **MORTISE TENON JOINTS**
Interlocking puzzle-like joints

STEEL FRAME | ■
Heavy-duty metal framework is forged to last.

..... ■ | **HEAVY-DUTY LIFT ACTUATOR AND SCISSOR MECHANISM**



Design also includes integrated battery backup and quiet and smooth lift system.

1. Pride® Power Recliner Collection products are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



INFINITY¹

FDA Class II Medical Device¹





When it's been a long day and you need to rejuvenate your body, mind and spirit, consider the power of true-infinite positioning. The Infinity collection by Pride® Power Lift Recliners is designed to be fully adjustable to your comfort needs. At the touch of a button on the sleek hand remote, the backrest and footrest move independently, so you get the support you need wherever you need it.

With pocketed coils that conform to your body, feel the stress of the day fade away as you stretch out fully with the footrest extension*. The Infinity collection is the ideal addition to your living room or family room. Available in four different fabric selections, the Infinity Recliner is built to last, year after year. As your life changes, one thing remains the same. The absolute comfort and support of an Infinity Power Lift Recliner.



STANDARD FABRIC: BURBANK



PEWTER

TAUPE

NAVY

JAVA

*Footrest extension is standard on all Infinity models except the 525iPW.



HERITAGE¹

FDA Class II Medical Device¹





There is a chill in the morning air as a new day begins. Before taking on the bustle of your life, embrace the inviting warmth and coziness of the Heritage collection by Pride® Power Lift Recliners. Stretch out completely with the footrest extension* and enjoy the steaming aroma from your morning coffee. Total Comfort Seating ensures outstanding relaxation and comfort, as you feel the softness of luxurious Cloud 9 fabrics enveloping your body.

When you can't put the day off any longer, innovative First Up® technology on the Heritage collection lifts you quickly and smoothly, helping you to get up and get moving. With a quiet and smooth lift system that offers durable and long-lasting performance, the Heritage collection fulfills your every need. For the recliner that perfectly complements every aspect of your life, choose the Heritage collection by Pride.



STANDARD FABRIC: CLOUD 9



BLACK CHERRY PACIFIC STONE WALNUT

*Footrest extension is not available on the LC-358S, LC-358XL & LC-358XXL.





LC-105

When it comes to the essentials of everyday life, never sacrifice your comfort. Experience a simply elegant solution with the Essential LC-105 Power Lift Recliner. Fill the empty spot in your house or apartment with the chair that complements your décor perfectly. With beautiful micro-suede fabrics in soft, dark colors, the Essential Lift Chair from Pride® Power Lift Recliners is the accent piece you've been waiting for!

The Essential LC-105 has an easy-to-use hand remote, so finding a comfortable position is as simple as pressing a button! Sinuous Spring Seating gives you endless support. Rolled, pleated armrests add a touch of elegance. With dual pockets that are easily accessible, the essentials of your life are always at your fingertips. Discover the elegant simplicity of the Essential LC-105 Recliner.



STANDARD FABRIC: MICRO-SUEDE



OPTIONAL FABRIC:
LEXIS URETHANE®

COCOA

SANDAL

SKY

NEW
CHESTNUT

ESSENTIAL¹



LC-250

Although the demands of life are relentless, take solace in the comfort of a Pride® Power Lift Recliner. The Essential LC-250 Recliner gives you the boost you need to take on the day. Total Comfort Seating provides complete support for ultimate relaxation. The two-button hand remote is simple to use. Adjust your chair effortlessly to find your perfect position every time.

With luxurious Cloud 9 fabrics that are soft to the touch, you'll find the perfect match to suit your home décor. Enjoy the plentiful storage of dual pockets that are easy to reach. Experience quiet and smooth lift performance, day after day, year after year. It's time to invest in your comfort with the Essential LC-250 Recliner.



STANDARD FABRIC: CLOUD 9




BLACK CHERRY



PACIFIC



STONE



WALNUT



OPTIONAL FABRIC: LEXIS STA-KLEEN®



CHESTNUT

ESSENTIAL¹

OPTIONAL FABRICS

ULTRAFABRICS® ULTRALEATHER®



BUFF FUDGE GARNET PECAN CHARCOAL

Ultrafabrics® Ultraleather® on Pride® Power Lift Recliners offer a richly detailed grain pattern expressive of premium grade leather. Ultrafabrics polyurethane products are manufactured utilizing Tukumi™ technology, a 4-layer laboratory process to create a luxurious, durable and comfortable polyurethane fabric.

FEATURES AND BENEFITS

- Rich feel and natural grain pattern expressive of soft luxury-grade leather
- Clinically proven more resistant to temperature changes than vinyl and genuine leather for comfort
- Antimicrobial surface cleans effortlessly with soap and water

CRYPTON® SUPER FABRIC




ESPRESSO COOL GREY RED SAND

Crypton® Super Fabrics from Pride® Power Lift Recliners are not only luxurious, durable and contemporary, but also resist spills, stains, odors and bacteria. Nothing gets past Crypton's unique patented integrated fiber process, protecting your upholstery inside and out.

FEATURES AND BENEFITS

- 100% water resistant – spills just wipe away
- Only disinfectable fabrics in the U.S. approved by the EPA
- Permanent technology that will not come off the fabrics
- 5-year warranty on moisture barrier

LEXIS STA-KLEEN® PERFORMANCE FABRIC



CHESTNUT MUSHROOM BLACK BURGUNDY

Lexis Sta-Kleen® Performance Fabric is a revolution in faux leather upholstery — the first stain resistant urethane fabric. Sta-Kleen offers permanent protection, created with a proprietary bonding process inherent in the manufacturing process of the product. It is not a topical application; it will last the lifetime of the upholstery.

FEATURES AND BENEFITS

- Look and feel of real leather, but with care-free durability
- Exceptional strength, abrasion resistance, and UV stability
- "Stain shield" allows spills to wipe off with a cloth, and it's anti-bacterial

QUICK REFERENCE FABRIC GUIDE

| MODEL | BURBANK | CLOUD 9 | MICRO-SUEDE | LEXIS STA-KLEEN | CRYPTON ARIA | ULTRALEATHER | LEXIS URETHANE |
|-----------|----------|--------------------------------|-------------|-------------------------------------|--------------|--------------|----------------|
| LC-525iPW | Standard | N/A | N/A | Optional | Optional | Optional | N/A |
| LC-525iS | Standard | N/A | N/A | Optional | Optional | Optional | N/A |
| LC-525iM | Standard | N/A | N/A | Optional | Optional | Optional | N/A |
| LC-525iL | Standard | N/A | N/A | Optional | Optional | Optional | N/A |
| LC-358PW | N/A | Standard | N/A | Optional | Optional | Optional | N/A |
| LC-358S | N/A | Standard | N/A | Optional | Optional | Optional | N/A |
| LC-358M | N/A | Standard | N/A | Optional | Optional | Optional | N/A |
| LC-358L | N/A | Standard | N/A | Optional | Optional | Optional | N/A |
| LC-358XL | N/A | Standard | N/A | Optional | Optional | Optional | N/A |
| LC-358XXL | N/A | Standard (only Stone & Walnut) | N/A | Optional (only Chestnut & Mushroom) | N/A | N/A | N/A |
| LC-105 | N/A | N/A | Standard | N/A | N/A | N/A | Optional |
| LC-250 | N/A | Standard | N/A | Optional (only Chestnut) | N/A | N/A | N/A |

ACCESSORIES



HEAD & LUMBAR PILLOW

(SOLD SEPARATELY)⁵

Adds additional support and comfort (in the same fabric the chair was ordered in) in the head or lumbar area.



BOLSTER PILLOW

(AVAILABLE ON MOST PRIDE® POWER LIFT RECLINERS)⁵

Fabrics to match your chair. Round, firm pillow with zipper to adjust the fill. Adds additional customized support where needed.



NON-SKID LEG LEVELERS

Non-skid leg levelers keep your Power Lift Recliner from sliding on a vast array of flooring.



DELUXE HEAT AND MASSAGE/ HEAT AND MASSAGE⁴

Deluxe heat and massage features four massagers and two heat pads. Regular heat and massage offers two massage units and one heat pad.



FIRST UP® TECHNOLOGY¹

Patent-pending technology offers up to a 50% increase in lift speed vs conventional lift chairs.



CLEANING KIT

This kit can be used with any Power Lift Recliner fabric (excluding vinyl/leather) in order to remove almost any stain or odor.



LITHIUM BATTERY BACKUP

(NOT AVAILABLE ON LC-105 & LC-250)

Next time the electricity goes out, be prepared with the new optional rechargeable lithium backup battery.

| COLLECTION | INFINITY ¹ | INFINITY ¹ | INFINITY ¹ | INFINITY ¹ |
|---------------------|--|---|---|---|
| MODEL | LC-525iPW | LC-525iS | LC-525iM | LC-525iL |
| POSITION | True-Infinite Position Zero Gravity 5'3" & Below | True-Infinite Position Zero Gravity 5'3" & Below | True-Infinite Position Zero Gravity 5'4" - 5'9" | True-Infinite Position Zero Gravity 5'9" - 6' |
| SIZE | Petite Wide | Small | Medium | Large |
| WEIGHT CAPACITY | 375 lbs. | 375 lbs. | 375 lbs. | 375 lbs. |
| SEAT-TO-FLOOR | 19" | 19" | 20" | 20" |
| SEAT DEPTH | 19" | 19" | 20" | 21" |
| SEAT WIDTH | 24" | 19.5" | 20" | 21.5" |
| TOP OF BACK TO SEAT | 28" | 25" | 28" | 29" |
| FULLY PADDED CHAISE | Yes | Yes | Yes | Yes |
| BACK STYLE | Blown Fiber Split Back | Blown Fiber Split Back | Blown Fiber Split Back | Blown Fiber Split Back |
| ADDITIONAL OPTIONS | Deluxe Heat & Massage ⁴ , Lithium Battery Backup | Deluxe Heat & Massage ⁴ , Lithium Battery Backup | Deluxe Heat & Massage ⁴ , Lithium Battery Backup | Deluxe Heat & Massage ⁴ , Lithium Battery Backup |
| FOOTREST EXTENSION | N/A | Standard | Standard | Standard |
| SEATING TYPE | Pocketed Coil Springs | Pocketed Coil Springs | Pocketed Coil Springs | Pocketed Coil Springs |
| DISTANCE FROM WALL | 19" | 18.5" | 19" | 21.5" |
| OVERALL WIDTH | 38" | 34" | 35" | 37" |
| CHAIR WEIGHT | 130 lbs. | 131 lbs. | 134 lbs. | 137 lbs. |
| WARRANTY | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated |

*Not all units achieve all positions. Please refer to specific product.

1. Pride® Power Recliner Collection products are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

2. Available on Heritage Collection except LC-358XL & XXL | 3. Available on LC-358 S, PW, M & L | 4. USB hand control not available with heat & massage options. | 5. Not available on all models.

| HERITAGE ¹ | HERITAGE ¹ | HERITAGE ¹ | HERITAGE ¹ | HERITAGE ¹ | HERITAGE ¹ | ESSENTIAL ¹ | ESSENTIAL ¹ |
|---|---|---|---|--|----------------------------|----------------------------|---------------------------|
| LC-358PW | LC-358S | LC-358M | LC-358L | LC-358XL | LC-358XXL | LC-105 | LC-250 |
| 3-Position 5'3" & Below | 3-Position 5'3" & Below | 3-Position 5'4" - 5'9" | 3-Position 5'9" - 6' | 3-Position 5'8" - 6'1" | 2-Position 5'8" - 6'1" | 3-Position 5'3" & Below | 3-Position 5'4" - 5'9" |
| Petite Wide | Small | Medium | Large | Heavy Duty | Very Heavy Duty | Small | Medium |
| 375 lbs. | 375 lbs. | 375 lbs. | 375 lbs. | 500 lbs. | 600 lbs. | 325 lbs. | 375 lbs. |
| 19" | 19" | 20" | 21" | 20" | 20" | 18" | 20" |
| 19" | 19" | 20" | 22" | 22" | 22" | 19" | 20" |
| 24" | 18" | 20" | 22" | 26" | 30" | 20" | 22" |
| 23.5" | 24" | 24.5" | 25.5" | 29" | 29" | 28" | 26" |
| Yes | Yes | Yes | Yes | Yes (removable) | Yes (removable) | Yes | Yes |
| Blown Fiber Button Back | Blown Fiber Button Back | Blown Fiber Button Back | Blown Fiber Button Back | Blown Fiber Button Back | Blown Fiber Button Back | Sewn Pillow Back | Blown Fiber Split Back |
| Heat & Massage ⁴ , First Up [®] , Lithium Battery Backup | Heat & Massage ⁴ , First Up [®] , Lithium Battery Backup | Heat & Massage ⁴ , First Up [®] , Lithium Battery Backup | Heat & Massage ⁴ , First Up [®] , Lithium Battery Backup | Heat & Massage ⁴ , Lithium Battery Backup | Lithium Battery Backup | No | No |
| Standard | N/A | Standard | Standard | N/A | N/A | N/A | N/A |
| Total Comfort | Total Comfort | Total Comfort | Total Comfort | Sinuuous Spring | Sinuuous Spring | Sinuuous Spring | Total Comfort |
| 18" | 18" | 17" | 19" | 19" | 14" | 17" | 20" |
| 40" | 35" | 35" | 37" | 39.5" | 43" | 33" | 34" |
| 130 lbs. | 118 lbs. | 117 lbs. | 120 lbs. | 159 lbs. | 168 lbs. | 96 lbs. | 118 lbs. |
| 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 1-Yr Limited | 2-Yr. Limited |

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; We reserve the right to alter specifications without prior notice.



Look for the Quick Ship Indicator for the shortest lead times available.

PRIDE[®]
POWER
RECLINER
COLLECTION¹

(US) 800-800-8586 | (CAN) 888-570-1113 | pridemobility.com |    

