

TRANSFORMER 2

FOLDING SCOOTER

Ultra Lightweight

Remote Operated Folding Scooter

Meets FAA Guidelines for Airline Travel



EMobility

Enhance Your Life



LED Headlight

Updated Control Pod

Wireless Remote

Lightweight Lithium Ion Battery

Specifications	Transformer S3026
Overall Length	38"
Overall Width	22" w/ Armrests
Tires	Front – 6x1.5" Rear – 7x2.36"
Maximum Speed	3.75mph
Range	Up to 13.5 miles
Turning Circle	55"
Weight (w/ battery)	55lbs
Color	Red, Blue, Yellow, Black
Brake System	Intelligent Regenerative Electromagnetic Braking
Drive System	Rear Wheel Drive
Ground Clearance	3 1/4 "
Weight Capacity	Up to 300lbs
Seat Width	17"
Seat Height (from ground)	22"
Battery	24V10Ah Lithium Ion
Motor	120W
Battery Charger	DC24V/2A
Folded Dimensions	18" x 22" x 25"



Enhance Your Life
 For more Information...
www.EnhanceMobility.com
 1-800-296-4151
Sales@EnhanceMobility.com